

**AED Plus®**  
Accessories and Training Products

**ZOLL®**



**ACCESSORIES FOR EVERY RESCUE**



### 7-inch Standard Wall Cabinet

7-inch alarmed wall cabinet holds AED Plus® (without carry case) and one spare set of electrodes  
**Part# 8000-0817**

## WALL CABINETRY

Solutions to protect your AED and improve visibility



### Semirecessed Wall Cabinet

- Alarmed AED wall cabinet holds AED Plus (without carry case) and one spare set of electrodes
- Meets ADA (American Disability Act) guidelines
- Professional installation recommended

**Part# 8000-0814**



### Fully Recessed Wall Cabinet

- Alarmed AED wall cabinet holds AED Plus (without carry case) and one spare set of electrodes
- Meets building code guidelines that prohibit cabinet from protruding more than 2 inches from wall
- Meets ADA (American Disability Act) guidelines
- Professional installation recommended

**Part# 8000-0811**



### 9-inch Standard Wall Cabinet

- 9-inch alarmed wall cabinet accommodates AED Plus in carry case with one spare set of electrodes

**Part# 8000-0855**

# ELECTRODES, BATTERIES & ACCESSORY KIT

Everything you need to keep your AED Plus ready to use



## CPR-D-padz® Electrodes

- One-piece design ensures fast and accurate placement
- For use on adult patients and children over the age of 8 or > 55 pounds (25 kilograms)
- Built into the electrodes, Real CPR Help® technology senses the rate and depth of chest compressions
- Five-year shelf life

Part # 8900-0800-01



## Pedi-padz® II Pediatric Electrodes

- ZOLL's intelligent pediatric capability is built into the AED Plus. The AED can sense when Pedi-padz II are being used and will switch the AED Plus from adult to pediatric mode, modifying the device so it will provide an appropriate pediatric ECG analysis algorithm and therapy energy level.
- For use with patients under the age of 8 years or < 55 pounds (25 kilograms)
- Two-year shelf life

*Caution: Pediatric electrodes should never be used with an adult patient*

Part # 8900-0810-01



## Accessory Kit

- The AED Plus accessory kit includes scissors, disposable razor, face barrier/shield, antimicrobial wipe, and dry towel.

*Note: This kit is included and attached to every set of ZOLL CPR-D-padz® electrodes.*

Part # 8900-0807-01



## Type 123 Lithium Batteries (sleeve of 10)

- Nonproprietary off-the-shelf Duracell 123 lithium batteries
- 10-year shelf life
- 5-year standby life once installed in the AED Plus

*Caution: Only use batteries from a ZOLL-approved manufacturer (refer to User Manual).*

Part # 8000-0807-01

# MOUNTING BRACKET & SIGNAGE

Ensure AED visibility



## Mounting Bracket

- Bracket is used to mount AED Plus on wall without carry case

Part# 8000-0809-01



## AED Wall Decal

- 4" x 4" AED decal suitable for indoor or outdoor use

Part# 8000-0849-01



## AED Sign Package

- 3D wall sign and 8.5" x 11" wall-mounted sign to improve AED visibility

Part# 8000-0825

# CARRY CASES & BAGS

Options for storage and transport



## Large Pelican Case

- Protection for AED Plus in more demanding environments
- Large enough to store AED Plus, CPR-D-padz, Pedi-padz II, and sleeve of 10 batteries
- Dimensions: 19" x 15.4" x 7.6" (48.3 cm x 39.1 cm x 19.3 cm)
- IP67 certified
- Military Standard C-4150J

Part# 8000-0837-01



## Small Pelican Case

- Protection for AED Plus in more demanding environments
- Stores AED Plus unit only
- Dimensions: 16" x 13" x 6.9" (40.6 cm x 33 cm x 17.5 cm)

Part# 8000-0836-01



## Soft Carry Bag

- Provides easy access and quick retrieval of the AED Plus and spare accessories
- Pockets for storing extra electrodes and batteries

Part# 8000-0802-01



# AED DEVICE MANAGEMENT & UPGRADE TOOLS

## RescueNet® Code Review Software

- Software providing ECG, vital sign data, and rate and depth of compressions for complete rescue review
- Can be downloaded at no cost directly from ZOLL website

**Part# 8000-0608-01**

## Administrative Software (ZAS)

- Software to reconfigure settings on ZOLL AED Plus
- Requires use of infrared IrDA adapter (sold separately)
- Software can be downloaded at no cost directly from ZOLL website

**Part# 8000-0311-01**

## AHA Guidelines Upgrade Kit

- Kit upgrading ZOLL devices manufactured prior to 2010 to current AHA Guidelines
- Includes Administrative Software (ZAS) CD, instructions, and labels to upgrade AED Plus faceplate
- Requires use of IrDA adapter (sold separately)

**Part# 7771-000011-01**

## USB-to-IrDA Adapter Kit

- USB-to-IrDA adapter cable and driver software
- Required for AED Plus configuration and post-rescue event data review

**Part# 8000-0815**

## Mini USB-to-IrDA Adapter Kit

- Mini USB-to-IrDA adapter cable and driver software
- Required for AED Plus configuration and post-rescue event data review

**Part# 8000-1200**

## RS232-to-IrDA Adapter Kit

- RS232-to-IrDA adapter cable and driver software
- Required for AED Plus configuration and post-rescue event data review

**Part# 8000-0816**

# AED PLUS CLINICAL UNIT TRAINING DEVICES & ACCESSORIES



## AED Plus Demo Kit

- Kit to demonstrate Real CPR Help® with clinical AED Plus
- Includes CPR-D Demo electrodes and manikin in large carry bag
- Clinical AED Plus and simulator not included

Part # 8000-0834-01



## CPR-D Demo Electrodes

- Demonstration electrode pads for use with clinical AED Plus
- Requires simulator

Part # 8900-5007



## Demo Manikin

- Adult manikin for AED/CPR training
- Use with AED Plus clinical unit and AED Plus Trainer2
- Training pads and simulator sold separately

Part # 8000-0835-01



## AED Plus Simulator

- Simulates heart rhythm for AED Plus demonstration using clinical unit

Part # 8000-0819-01



## Travel Trainer

- One-piece, compact, easy-to-carry demonstration tool for clinical AED Plus
- Simulates VF and normal heart rhythms

Part # 8008-0006-01

# AED PLUS TRAINER2 & ACCESSORIES

Everything you need to simulate a rescue with your AED Plus Trainer2



## AED Plus Trainer2

- Allows the student to experience what it's like to help save a life with the AED Plus. Choose from four different training scenarios or have complete manual control of the entire rescue. A wireless remote lets you control multiple Trainer II devices.
- Ships with Operator's Manual, batteries, wireless remote, and adult training electrodes
- Replacement adult training electrodes, pediatric training electrodes, and replacement adhesive gel can be purchased separately as needed

Part # 8008-0050-01 (semiautomatic)

Part# 8008-000052-01 (automatic)



## Pedi-padz II Training Electrodes

- Pediatric electrodes for use with the Trainer2

Part # 8900-000861-01  
(pack of 6)



## Replacement Trainer II Wireless Remote

- Wireless remote for use with the Trainer2

Part # 8008-0007



## Replacement Adhesive Gel for CPR-D-padz Training Electrodes

- Gel must be applied to the reusable pad before use.

Part # 8900-0803-01 (pack of 5)



## AED Plus Trainer2 Soft Carry Case

- Clearly identifies the device as a training unit and provides a convenient way to transport the Trainer2. Includes a pocket for accessories.

Part # 8000-0375-01



## Demo Manikin

- Adult manikin for AED/CPR training use with AED Plus clinical unit and AED Plus Trainer2
- Training pads and simulator sold separately

Part # 8000-0835-01



## CPR-D-padz Training Electrodes

- Adult electrodes for use with the Trainer2. Replacement gel can be purchased separately.

Part # 8900-0804-01



# AED PROGRAM MANAGEMENT

Stay prepared to save a life

## PlusTrac Professional AED Program Management for a single AED

- Medical Authorization/Rx
- Online management software access
- Protocol and EMS templates
- AED Law library and support knowledge base
- Email reminders and alert notifications
- Complete customer support/implementation
- Post-event services and loaner program

**1 Year**

**Part# 8000-1110-01**

**5 Year**

**Part# 8000-1111-01**

## PlusTrac Essentials AED Program Management for a single AED

- Medical Authorization/Rx
- Online management software access
- Protocol and EMS templates
- AED Law library and support knowledge base
- Email reminders and alert notifications

**1 Year**

**Part# 8000-1116-01**

**5 Year**

**Part# 8000-1117-01**

For general product-related questions or assistance with where to buy ZOLL accessories, contact Customer Support at 800 480-1277. For technical questions about your ZOLL AED Plus or accessories, contact Technical Support at 978 421-9460.

## CPR SAVERS & FIRST AID SUPPLY

7904 E. Chaparral Rd. | Scottsdale, AZ 85250 | 800-480-1277 | [www.cpr-savers.com](http://www.cpr-savers.com)

**ZOLL**®