

Meds Pack 32375

> Dimensions: 13"L x 3"W x 5.5"H

Iron Duck's Ultra Meds Pack is a compact drug module allowing separate storage, easy restocking and quick access to delicate ampoules and vials, whether boxed or unboxed.

The exterior shell is constructed of 1000D Nylon and features a web quick-grab handle. The front features a full-length clear vinyl document window. The back features a full length strip of loop to secure the module into any bag or pack with hook strips on the liner.

The case contains foam inserts with die cut grids that allow for a custom fit of medication containers. Clear vinyl windows close over medications to hold them in place. Each foam insert is secured in place by hook and loop to facilitate easy replacement.

