



## Smart Pack Airway Backpack 32410

Overall Dimensions : 19" h x 15" w x 6.5" d

Constructed from 1000 denier nylon the Smart Pack Airway Backpack will help you carry your airway supplies no matter where your day takes you.

Open this durable pack to find two inch nylon straps to secure either a "D" or Jumbo "D" oxygen tank. On each side of the tank you will find a pouch for different needs. Above the tank straps, you will find a large utility pouch featuring a clear vinyl window to view your contents with ease. Below the tank you will find a quick draw mesh pouch. Each of these pouches are secured with a hook and loop system and can be removed with ease.

Inside the front panel you will find a variety of pockets and elastic loops to help keep you organized. To help keep everything clean there is a zippered vinyl cover that protects your crucial supplies. The face of the backpack has a large pouch to slide in papers or other flat supplies. Throw the Smart Pack Airway Backpack over your shoulder with the single shoulder strap or carry it under your arm.