



## IDX Expandable Drug Module 34510

- > Standard Nylon Colors: Red
- > Dimensions: 9"L x 6.5"W x 3.5"D (6"D fully expanded)

Iron Duck's IDX Expandable Drug Module supports easy restocking and fast deployment via notebook style binding. The Drug Module can hold up to 8 removable & customizable pages with multiple pockets and elastic tacking for supplies.

The exterior shell is made of 1000 denier nylon with dual, reverse-coil YKK zippers. On the front there is a clear vinyl window for easy identification of medications and restocks dates. The unique, expandable zipper increases storage capacity by 175.5 cubic inches. Pull open the heavy duty zippers to take a look at the multitude of organization tools inside the main compartment. The left side of the main compartment offers a flat zippered pocket layered under two stacked flat pockets for alcohol preps and dressings. A tape dispenser completes the left side. On the right side you will find a clear vinyl gusseted pocket beneath an elastic pocket. A row of elastic tacking holds used syringes or other tools. In the middle, you will find four pages that you can customized to your needs.

The reverse side of page has a flat zippered pocket, clear window for documents and a row of elastic tacking for syringes. Each page contains the following details:

- \* 34510-P1 – 2 rows of 6 clear pockets with a flat zippered pocket on the opposite side with elastic tool holders
- \* 34510-P2 – 2 rows of 5 clear pockets with a flat zippered pocket on the opposite side with elastic tool holders
- \* 34510-P3 – 1 row of 5 clear pockets with a flat zippered pocket on the opposite side with elastic tool holders
- \* 34510-P4 – 2 clear pockets with "pick & pluck" foam blocks with a flat zippered pocket on the opposite side with elastic tool holders.

Carry your midwife Breathsaver with the heavy duty wrap around handles or with the agronomical, adjustable shoulder strap.