



## IDX Intubation Module 34520

- > Dimensions: 13"L x 8"W x 2.5"H
- > Dimensions open: 13"L x 27.5"W x 2.5"H)

Iron Duck's tri-fold Airway Module organizes a complete set of intubation supplies including laryngoscope handles and blades, airways, ET tubes and holders and accessories. The exterior shell is made of 1000D Nylon and features a quick-release plastic buckle closure and a clear vinyl card-sized window. A clear vinyl document panel on the back side completes the exterior design.

The inside left panel has a flat zippered pocket beneath two rows of elastic loops. The center panel holds a full set of airways in elastic loops. Clear gusseted zipper pocket holds ET tube holders and laryngoscope bulbs. A clear vinyl flat pocket is layered over the gusset for topical anesthetic and packets of lubricating jelly. A row of elastic loop holds two sizes of intubation forceps, multiple sized syringes and other small tools.

A tape dispenser completes the center panel. The inside right panel is configured with two clear vinyl pockets to hold ET Tubes, Adult and Pedi styles and bite sticks.

