


CAPACITY:
Pediatric Patients Aged 0-5 Years

DIMENSIONS:
Dimensions: 48" x 12" x 1.75"

X-RAY TRANSLUCENT:
100% X-Ray, MRI and CT Scan compatible.

DESIGN:
Patented Dual Plane Headrop System *(Fig.1)
Rotomolded one piece seamless HDPE shell over non-toxic Eco-Friendly polyurethane foam*

STANDARD COLORS:
Yellow, Teal



Patented Dual Plane Headrop System.

*(Fig.1)



*Pedi~Air~Align Shown with Straps.

Pedi Air Align

Product#: 35840 (Complete)
Product#: 35841 (Board Only)

The Pedi-Air-Align's patented dual plane head drop system allows for a patient to be completely immobilized while managing their airway and establishing neutral spinal alignment simultaneously. The board accommodates the larger occipital area of a child's head, allowing for neutral spinal alignment and an open, maintainable airway. The Pedi-Air-Align is the only pediatric immobilization board that is rotomolded from high-density polyethylene (HDPE) creating a one piece seamless product. This spine board is completely impervious and as a result is easy to clean and disinfect.

The Pedi-Air-Align Complete System includes:

1 board, 1 shoulder restraint, 1 chest strap, 4 body straps, 2 wrist restraints, 6 disposable head blocks, 6 disposable head block straps, and 1 carry case.

