

WEIGHT CAPACITY:

500 lbs.

WEIGHT:

16lbs.

DIMENSIONS:

72" x 16" x 1.75"

FEATURES:

250 lbs. Float Capacity
100% X-Ray Translucent,
MRI and CT Scan compatible

DESIGN:

1 Rotomolded one piece
seamless HDPE shell over non-toxic
Eco-Friendly polyurethane foam*
Fourteen over sized hand holds
• Available with 12 speed clip pins

RELATED PRODUCTS:

Straps & Strapping Systems
Sof Loc Head Immobilizer

STANDARD COLORS:

Yellow, Blue, Orange, Red
Lime Green, Olive Drab, Black



Base Board

Product#: 35850

With Pins: 35850-P

Iron Duck's BASE Board is the most rigid backboard in its weight class... bar none. Recent tests prove that this simple, effective and economical immobilization solution provides the lowest deflection in the 400 to 500 lb weight class! Flat patient surface won't compromise CPR efforts or aggravate cervical spine injuries (as concave boards may) while tapered design and convex bottom facilitate extrication, log-rolling and easy lifting.

The board is fabricated with a tough, completely impervious rotomolded HDPE shell over non-toxic eco-friendly polyurethane foam. Subtle surface texture prevents slippage but is easy to clean and disinfect. Fourteen oversized hand-holds each measuring 2" x 6" will accommodate gloved hands. 12 optional offset composite pins can be located in the lower third of the handholds, keeping pins out of the way when lifting patients.

This highly rigid, 100% X-Ray translucent spine board is easy to maintain and store, making it the only choice for economical daily use and mass-casualty.