

WEIGHT CAPACITY:

1000lbs

DIMENSIONS:

72" x 16" x 1.75"

WEIGHT:

14 lbs.

FEATURES:

100% X-Ray, MRI and CT Scan compatible
250 lbs. Float Capacity

DESIGN:

Rotomolded one piece seamless
Fourteen over sized
raised hand holds
Available with 12 speed clip pins

RELATED PRODUCTS:

Straps & Strapping Systems
Sof Loc Head Immobilizer

STANDARD COLORS:

Yellow, Blue, Orange, Red
Lime Green, Olive Drab, Black



Ultra Loc

Product#: 35900

With Pins: 35900-P

Extremely thin at 1.75" and lightweight at 14 lbs.

The Ultra-Loc by Iron Duck exhibits less deflection at applied weights up to 1,000 lbs. which is the top weight capacity in the market. The open leg triangle integrates with full-length, molded-in runners to offer the rigidity that concave boards lack while keeping the hand-holds off the ground for ease in lifting.

The flat patient surface won't compromise CPR efforts or aggravate cervical spine injuries while tapered design and molded-in runners facilitate extrication, logrolling and easy lifting. Fourteen oversized hand-holds each measuring 2" x 6" will accommodate gloved hands. Fabricated with a tough, completely impervious, rotomolded HDPE shell. The subtle surface texture prevents slippage but is easy to clean and disinfect.