



Oxygen Bag 1-Pocket 36002D-PK

Includes shoulder strap

Dimensions: 21"L x 11.5"W x 8.5"H ~ 1,989 cu in

This clamshell design oxygen bag with 1 full length exterior pocket provides added storage and is ideal for protection and fast deployment of a D or Jumbo D tank, regulator and accessories. Fully lined with high impact foam the removable liner makes cleaning easy. The main compartment is accessed through heavy duty zippers that run the length of the bag and down both ends.



A padded regulator protector and adjustable tie-downs hold the tank securely in place during transport. The exterior full length pocket contains multiple elastic holders for organizing supplies. Several interior pockets provide storage for masks, tubes and accessories.



An extra layer of UP fabric is added to the bottom of this bag for extra durability and ease of cleaning. Heavy duty handles wrap around the entire bag for improved weight distribution.