


Soft Stretcher

Product#: 38032



Alignment Strap



The Soft Stretcher allows you to have a high level of flexibility with your response need. Use this sheet to move patients, as means of transfer to other patient transport equipment or along with your backboard of choice.

The contrasting center stripe will help you center your patients without having to shift them around. Move your patient with ease by using the eight ergonomically designed rugged handles made out of heavy-duty webbing that will support patients of all sizes.

- Heavy-duty vinyl construction maximizes durability under harsh conditions
- Flexible material facilitates easy transfer for any size patient
- 8 rigid heavy-duty lift handles provide superior versatility and operator comfort
- The transfer-flat rolls into compact storage size with convenient carry handle
- Maximum Weight Capacity 1600 lbs



Alignment Strap